METAL-CERAMIC-ALLOY

ADORBOND CC PLUS

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METAL-CERAMIC-ALLOY

ADORBOND CC PLUS

ADORBOND CC PLUS is a dental metal-ceramic alloy based on Cobalt. ADORBOND CC Plus is free of nickel, beryllium, indium and gallium and fulfills the of the standards EN ISO 22674 type 4, for dental alloys. Through the low Vickers hardness of 320 HV10 result the alloy easy to drill, mill and polish and is therefore particularly well suited for primary situations.

ADORBOND CC PLUS has CE 0197 certification.

Composition		Technical Data (estimated values)	
(Mass-%)		Density (g/cm ³)	8,2
Cobalt	63,6	Vickers hardness (HV 10)	320
Chromium	24,8	Thermal expansion coefficient	
Molybdenium	5,0	25 - 500 °C (10 ⁻⁶ K ⁻¹)	14,3
Tungsten	5,5	20 - 600 °C (10 ⁻⁶ K ⁻¹)	14,5
Silicium	1,1	Melting range (°C)	1.360 – 1.420
		Casting temperature (°C)	ca. 1.500
		Yield strength (MPa)	460
		Modulus of elasticity (N/mm ²)	ca. 200.000
		Ductile yield A ₅ (%)	12

Waxing-up:

Wax-up with crown- and bridge- wax as usual. Prevent thickness of material lower than 0.35 mm. Lead wax sprues indirectly. For wax sprues use round wax wires with Ø 2.0-2.5 mm for single crowns and 2.5-3.0 mm for bridges. For frames with more than 4 teeth use a distribution funnel with Ø 3.5-4.0 mm, for massive pontics to Ø 5mm use.

Melting and Casting:

Suitable are phosphate bonded investments for crown and bridge work. Preheat the investment to about 850 to 900 °C. Hold temperature for about 30 minutes. Refer to manufacturer's instructions for use for the casting machines. For ADORBOND CC Plus use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after each use to avoid residues of slag. When melting by induction heating start casting as soon as the ingots have collapsed giving an uniform melt. For melting by flame heat the ingots and give a rotary motion by use of the flame. Start casting as soon as the bath begins to move. Allow the cylinder slow air cool down to the ambient temperature and deflask.

Preparing of Metal-Ceramic Restorations:

Use common available ceramics for base metal alloys like Vita (Omega, VM13), Ivoclar, DeTrey oder Ducera having a coefficient to match that of the alloy, being suitable for firing temperatures between 900°C and 980 °C.

Firing of the ceramic:

- Sand blast** the surface by use of a **pencil-blaster**. Blast all of the surface on which the porcelain is to be applied with aluminium oxide with a grain size of min **100** or preferably **250 µm**.
- Ultrasonically clean** the frame for **5 to 10 minutes** in distilled water or degrease with steam or with ethyl acetate.
- The oxide firing is optional, to be done at about **960 °C** for abt. **5 minutes**.
- Always remove the oxide layer after oxide firing by sand blasting with 120 or 250 µm of aluminium oxide and degrease again!**
- The opaque** is then applied on the surface by a first thin wash firing and a second evenly covering opaque layer.
Before firing always let **dry for 5-10 minutes** at abt. **600 °C!**
- Firing has to be done according to the ceramic manufacturer's instructions.
- After every firing step (dentine bake, build-up and glazing) **cooling phase until ca. 750 °C**.

Finishing:

After firing of the ceramic, polish the frame with suitable grinding and polishing instruments for base metal alloys up to high gloss.

Soldering and Welding:

Soldering before firing of the frame can be carried out with **ADOR CC cobalt-based solder** and high temperature flux. For welding with laser use **ADOR Co-Cr laser wire**.

Final Safety Notes

Metal dusts in principle are harmful. Use a dust extractor. Consider allergic hypersensitivities for technical personal and Patient to contents of the alloy. In case of suspected incompatibility with individual elements of this alloy, this should not be used.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The user is responsible for the correct processing of the alloy itself.