


REF-LINE

Look and feel natural



M1

M3

M4

M6

K3

K4

K6

N3

N4

N6

W3

W4

W6

H3

H4

H6



V1,V3

V3,V4,U3

V4,U3,U4

V4,V6,V8,U4

V4,V3,U3

V6,V4,U4

V8,V6,U6

V3,V4,U3

V4,V6,U3,U4

V6,V8,U4,U6

V3,V4,U3

V4,V6,U3,U4

V6,V8,U4,U6

V3,V4,U3

V4,V6,U3,U4

V6,V8,U4,U6



Q1

R2,R4,Q1,Q3,F2,F4

R2,R4,Q3,F2,F4

R4,R6,Q5,F4,F6

R4,F4,Q3,R2,F2

R4,F4,Q3,R2,F2

R6,F6,Q5,R4,F4

R2,R4,Q3,F2,F4

R2,R4,Q3,F2,F4

R4,R6,Q5,F4,F6

R2,R4,Q3,F2,F4

R2,R4,Q3,F2,F4

R4,R6,Q5,F4,F6

R2,R4,Q3,F2

R2,R4,Q3,F2,F4

R4,R6,Q5,F4,F6

SHADES

BL1 - BL3

A1 - A2 - A3 - A3,5 - A4



B1 - B2 - B3 - B4

C1 - C2 - C3 - C4

D2 - D3 - D4

PoliDent





M1

a 40,8
b 8,7
c 7,5



M3

a 43,2
b 9,7
c 8,2



M4

a 45,3
b 9,9
c 8,5



M6

a 49,1
b 10,5
c 9,1



K3

a 44,1
b 9,9
c 8,3



K4

a 46,1
b 10,3
c 8,7



K6

a 49,6
b 11,1
c 9,3



N3

a 43,6
b 9,9
c 8,3



N4

a 45,4
b 10,9
c 8,6



N6

a 48,8
b 11,4
c 9,3



W3

a 43,4
b 9,3
c 8,2



W4

a 45,3
b 9,7
c 8,5



W6

a 48,7
b 10,6
c 9,2



H3

a 43,0
b 10,7
c 7,8



H4

a 44,9
b 11,2
c 8,2



H6

a 48,3
b 12,0
c 8,8





V1

a 29,8
b 8,0
c 4,3

V3

a 31,9
b 8,3
c 4,7

V4

a 33,7
b 8,7
c 4,9

V6

a 35,8
b 9,2
c 5,3

V8

a 38,6
b 10,5
c 5,6

U3


a 34,1
b 8,4
c 5,1

U4

a 35,9
b 8,7
c 5,4

U6

a 38,3
b 9,4
c 5,7



a
b
c
d

R2

a 30,6
b 32,8
c 10,6
d 9,8

R4

a 33,2
b 35,4
c 11,4
d 10,7

R6

a 35,0
b 38,0
c 12,3
d 11,3

Q1


a 29,1
b 31,6
c 10,2
d 9,4

Q3

a 31,9
b 34,6
c 11,2
d 10,3

Q5

a 33,8
b 36,7
c 11,9
d 10,9



PoliDent



F2


a 30,4
b 32,4
c 10,2
d 9,5

F4

a 32,9
b 35,1
c 11,1
d 10,2

F6

a 35,4
b 37,7
c 11,9
d 11,0



R



F

abradirani zobje
lower cusps
abradirani zubi

R

Q

podaljšani vratovi
extended cervical on premolars
izduženi vrat



REF-LINE *Look and feel natural*