


			
3	31	I1	T1, S1
	32	I2, I3, L2	T1, T2, S1
	33	I1, L2, I3	T2, T3
	34	I3, I4	T2, T3, T4
	35	I2, I4, I5, L2	T3, T4
	36	I4, I5, L6	T3, T4
	37	I5, I6	T4, T6
	37A	I6	T4, T6
	38	I6, I8, I7	T4, T6, S5
	39	I7	T6
0	32I	I2, I3, L2	T1, T2, S1
	34I	I3, I4	T2, T3, T4
	37I	I6	T4, T6
	01	I1	T1, S1
	02	I1, I3, L2	T1, T2, T3, S1
	03	I3, I4	T3, T4, S1
4	04	I4, I5, L6	T3, S5
	05	I5, I6, L8	T4, S5
	05A	L8, I6	T4, S5, T6
	06	I6, I7, L8	S5, T6
	41	I1	T1, S1
	42	I3, I4, L2	T2, T3, S3
4	43	I3, I4, L2	T3, T4
	43A	I3, I4, L6	T3, S3, S5
	44	I5, I6, L6	T4, S5
	45	I5, I6	T4, S5
	46	L8	T4, S5, T6



31
38,6
8,9
7,2



32
41,5
9,5
7,9



33
42,0
10,5
8,2



34
43,5
9,4
8,0



35
44,1
10,2
8,4



36
44,9
10,4
8,5



37
46,5
10,8
8,7



37A
47,0
11,0
8,8



38
49,5
11,8
9,3



39
51,2
12,6
9,4

BARVE / BOJE / COLOURS / FARBE /
COLOR / COR / COLORE

BL3

A1 - A2 - A3 - A3,5 - A4

B1 - B2 - B3 - B4

C1 - C2 - C3 - C4

D2 - D3 - D4



YO - Y1 - Y2 - Y3 - Y4 - Y5 - Y6

B1 - B2 - B3 - B4 - B5 - B6

G1 - G2 - G3 - G4 - G5 - G6

PoliDent





01
40,3
8,2
7,5



02
42,3
8,4
8,0



03
44,2
9,8
8,2



04
46,6
9,8
8,7



05
47,3
10,0
8,7



05A
48,9
11,9
9,0



06
52,0
13,1
9,7



41
40,4
9,0
7,5



42
43,4
9,4
8,4



43
45,8
9,0
8,2



43A
46,8
9,4
8,6



44
46,9
9,4
9,0



45
48,5
10,0
8,9



46
51,4
10,3
9,2

CERVICAL WIDER SHAPES ALSO
SUITABLE FOR IMPLANTS




32I
43,5
9,4
8,0



34I
46,2
9,9
8,5



37I
49,9
10,7
9,1



I1

31,3
7,5
4,6



I2

31,4
9,1
4,7



I3

33,5
8,3
5,0



I4

34,0
9,1
5,0



I5

34,4
9,1
5,3



I6

37,5
8,5
5,6



I7

38,2
11,9
5,8



L2

34,7
8,2
5,1



L6

38,5
9,1
5,9



L8

40,2
9,3
6,2



T1

28,5
30,0
8,6
7,8



T2

28,7
31,2
8,3
8,2



T3

29,2
31,1
8,7
8,7



T4

30,5
33,2
9,0
8,7



T6

33,0
35,0
9,5
9,3



S1

27,2
28,2
7,9
7,6



S3

31,1
32,6
8,5
8,2



S5

31,6
33,2
8,9
8,8



primodent®
For beautiful smile.